



## VersaQuat™

VersaQuat™ is an EPA-registered, RTU (ready-to-use) one-step disinfectant and cleaner which utilizes well-established and accepted quaternary ammonium chemistry with acceptable residue and safety profiles. With broad-spectrum, best-in-class efficacy claims (30 kill claims), VersaQuat™ delivers consistent quality and intuitive usage. Recommended for use as a standard disinfectant as part of a disinfectant rotation program, VersaQuat™ is ideal for use in laboratories, hospitals, home healthcare settings, and veterinary clinics and can also be used in offices, schools, gyms, and homes. Quick, effective, and convenient, VersaQuat™ is not abrasive and is compatible with various hard, non-porous surfaces.



| FEATURES                                     | BENEFITS   |
|--|--|
| Kills 99.99% of bacteria, viruses, and fungi | • Kills common microbiological contaminations on surfaces                              |
| Consistent quality and intuitive usage       | • Can use in a variety of applications<br>• Ideal for large surface cleaning           |
| Utilizes safe and effective chemistry        | • No strong odor provides good user acceptability<br>• No special disposal is required |
| Bottles double bagged                        | • Facilitates transfer disinfection into cleanroom                                     |
| Not classed as corrosive                     | • Can be used safely on all hard, non-porous surfaces                                  |
| Ready to Use                                 | • No mixing required   |

[www.berkshire.com](http://www.berkshire.com)

Contact: Tel 1 800 242 7000 / 1 413 528 2602  
[info@berkshire.com](mailto:info@berkshire.com)

Rev. BC\_4020\_09-2022

Americas

Tel 1 413 528 2602

[usa@berkshire.com](mailto:usa@berkshire.com)

Europe

Tel 44 (0) 870 757 2877

[enquiries@berkshire.uk.com](mailto:enquiries@berkshire.uk.com)

SE Asia

Tel 65 6252 4313

[enquiries@berkshire.com.sg](mailto:enquiries@berkshire.com.sg)

Japan

Tel 81 3 5827 2380

[master@berkshire.co.jp](mailto:master@berkshire.co.jp)

| Bacteria   | Contact Time |
|--|--------------|
| <i>Pseudomonas aeruginosa</i> (ATCC 15442)   | 2 minutes    |
| <i>Staphylococcus aureus</i> (ATCC 6538)   | 2 minutes    |
| <i>Salmonella enterica</i> (ATCC 10708)  | 2 minutes    |
| Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA) (ATCC 33592)   | 2 minutes    |
| Vancomycin Resistant <i>Enterococcus faecium</i> (ATCC 51559)  | 2 minutes    |
| Multi-Drug Resistant (MDR) <i>Acinetobacter baumannii</i> Ceftazidime and Gentamicin Resistant (ATCC BAA-1605) | 2 minutes    |
| Extended Spectrum Beta-Lactamase (ESBL) <i>Escherichia coli</i> (ATCC BAA-196)                                 | 2 minutes    |
| Viruses  | Contact Time |
| Avian Influenza Virus (H5N1) NIBRG-14  | 30 Seconds   |
| Severe Acute Respiratory Syndrome Related Coronavirus 2 USA-WA1/2020 BEI NR-52281                              | 30 Seconds   |
| Human Rotavirus (ATCC VR-2018)   | 1 minute     |
| Norovirus (Norwalk Virus)  | 2 minutes    |
| Canine Parvovirus  | 2 minutes    |
| Rhinovirus Type 14 (ATCC VR-284)   | 2 minutes    |
| Mildew   | Contact Time |
| <i>Aspergillus niger</i> (ATCC 6275)   | 10 minutes   |

## Product Application

- Controlled environments and general institutional/industrial cleaning: Broad-spectrum efficacy (bacteria, viruses, fungi) and one-step cleaner and disinfectant.
- Sensitive surfaces and equipment: Good surface compatibility offers versatility and is not damaging.
- Healthcare and long-term care facilities: One step, ready to use effectively and efficiently clean and disinfect near-patient surfaces and equipment.
- Hospitality and school environments: No harsh chemicals, toxins, or added fragrances provide flexibility for use.

## Order Information

| Product    | Number  | Size (Net Contents) | Style                    | Bottles/Case |
|------------|---------|---------------------|--------------------------|--------------|
| VersaQuat™ | VQUAT6Q | 32 Fl. Oz. (946mL)  | Bottle w/trigger sprayer | 6            |



## Recommended Wipes to use with this product:

| Product                         |
|---------------------------------|
| <a href="#">Durx® 670</a>       |
| <a href="#">Durx® 770</a>       |
| <a href="#">Polx® Nonwoven</a>  |
| <a href="#">ValuClean® Plus</a> |
| <a href="#">MicroSeal® - VP</a> |