

**Product Information** 

ISO Class 5 - 8

Cleanroom Class 100 – 100,000 EU Grade A/B/C/D

# Sterile SatPax® 550

Sterile Pre-wetted 100% Meltblown Polypropylene Cleanroom Wiper

Sterile SatPax® 550 wipers are composed of exceptionally clean100% meltblown polypropylene pre-wetted with 70% IPA and 30% DI water. This pre-wetted format provides a cost-effective and easy-to-use solution versus traditional bulk handling of solvents, maintenance of squirt bottles and inconsistent wetting and cleaning associated with wetting a dry wiper. Each package is gamma irradiated sterile to a 10-6 Sterility Assurance Level making the product ideal for cleaning critical sterile processing environments.



## Other Sterile Pre-wetted wipers

- Sterile SatPax® 670
- Sterile SatPax® 1200
- Sterile SatPax® MicroSeal® WFI

### **Key Attributes**

- 100% meltblown polypropylene
- · No chemical binders in base material
- Pre-wetted with consistent 70% IPA/ 30% DI Water
- Re-sealable solvent resistant packaging
- Wipers in c-folded configuration for single withdrawal
- Gamma irradiated and sterile validated to a 10<sup>-6</sup> Sterility Assurance Level per AAMI guidelines
- Certificate of Sterility with each case and available on the web by lot number 24/7. Details: Expiration Date, Radiation Dosage, Lot Information, Meets cGMP requirements for traceability

#### **Benefits**

- Very low extractables and ions
- · Very low particle and fiber counts
- Reduces alcohol usage and preparation / handling costs
- Reduces VOC emissions
- Increases cleaning efficiency
- Increases cleaning protocol consistency

#### **Environmental**

Registered under REACH

## **Applications**

- Designed for use in ISO Class 5 and higher sterile cleanroom environments
- Designed for use in wet cleaning of critical surfaces where control of flammable solvents and flammable solvent concentrations is required
- Final cleaning of surfaces or products prior to manufacturing or packaging

#### Saturation Levels

The amount of solution contained in each wiper can be varied according to customer requirements. Higher saturation levels apply more solution to the surface during cleaning.

#### **Validated Sterile**

ANSI/AAMI/ISO procedures have been followed for verification of irradiation dose with audits being performed to monitor bioburden levels and dose. The Sterility Assurance Level (SAL) or the probability of a viable organism being present on a product unit after sterilization is a SAL of 10<sup>-6</sup> (which is the chance of 1 non-sterile unit in 1,000,000).

## www.berkshire.com

Contact: Tel 1 800 242 7000 / 1 413 528 2602 info@berkshire.com

America	Tel 1 413 528 2602	info@berkshire.com
Europe	Tel + 44 1953 562800	enquiries@berkshire.uk.com
SE Asia	Tel 65 6252 4313	enquiries@berkshire.com.sg
Japan	Tel 81 3 4530 9883	master@berkshire.co.jp



## **Technical Data (In Dry State)**

Attribute		Units	Value	Test Method
Basis Weight		g/m²	35.4	TAPPI T-410
Caliper		μm	204	TAPPI T-411
Fibers	<u>≥</u> 100μm	fibers/cm <sup>2</sup>	0.48	IEST-RP-CC004.4 Sec 7.1.3 /Sec 7.2.2 modified
Particles	<u>≥</u> 0.5μm	x10³/cm²	4.8	IEST-RP-CC004.4 Sec 7.1.3 /Sec 7.2.1 modified
Sorbency	Capacity	mL/m²	196*	IEST-RP-CC004.4 Sec 9.1 /Sec 9.2 modified
	Efficiency	mL/g	5.5*	
	Rate	seconds	10*	
Non-Volatile Residue	DI Water	g/m²	0.0032	IEST-RP-CC004.4 Sec 8.1.2
	IPA	g/m²	0.015	
Ions	Na⁺	ppm	0.42	IEST-RP-CC004.4 Sec 8.2.2
	K <sup>+</sup>	ppm	0.20	
	Ca <sup>++</sup>	ppm	0.47	
	Mg <sup>++</sup>	ppm	0.038	
	CI-	ppm	0.80	

#### Notes:

- Technical data represented in this table are typical values at the time of publication. These should not be used as
  product specifications.
- Due to differences in test methods applied and equipment utilized by different wiper manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions
- Third party testing can be performed upon request
- \*Sorbency Capacity, Efficiency, and Rate are performed with 70% IPA

# **Order Information:**

Product	Number	Size	Shts/pk	Pks/cs	IPA/DI	VOC % by Weight	Style
Sterile SatPax® 550	SSP55000424	9x11" (23x28cm)	50	24	70/30	48%	C-fold
Sterile SatPax® 550	SSP55000524	7x11" (18x28cm)	50	24	70/30	48%	C-fold

# Other Berkshire Products



Wipers



**Glove Liners** 



**Mop Systems** 



**Documentation Systems** 



**Face Masks** 



Swabs