

Gamma Wipe® 67

55% Cellulose / 45% Polyester Sterile Nonwoven Cleanroom Wiper

Gamma Wipe®67 is a polyester / cellulose blend wiper recommended for ISO Class 5 and above environments composed of 45% polyester and 55% cellulose. Each package is gamma irradiated to a 10⁶ Sterility Assurance. This combination of properties provides Gamma Wipe®67 with the level of cleanliness, abrasion resistance and chemical compatibility required for cleaning general sterile processing areas where both economy and cleanliness are required



Key Attributes

- 55% cellulose / 45% polyester hydroentangled nonwoven
- No secondary binders
- Available in two sizes: 9x9" and 12x12"
- Gamma irradiated and sterile validated to a 10⁶ Sterility Assurance Level per AAMI guidelines
- Header style packaging with Package Integrity testing
- Lot number on each pack
- Certificate of Sterility with each case and available on the web by lot number 24/7. Details: Expiration Date, Radiation Dosage, Lot Information, Meets cGMP requirements for traceability

Benefits

- Excellent combination of the synthetic polyester strength and cleanliness with the absorbent characteristics of cellulose.
- Smooth, highly absorbent and durable with low extractable levels
- Economical
- Chemically compatible with common cleaning and disinfecting solutions including isopropyl alcohol, phenols, bleaches and quarternary ammonium compounds

Applications

- Designed for use in ISO Class 5 and higher sterile environments and USP <797> applications
- Designed for cleaning sterile surfaces and equipment during aseptic processing
- General wiping in component prep, compounding and wash areas
- General wipedowns in sterile gowning rooms
- Applying and removing cleaning and disinfecting solutions

Pre-Wetted Option

The same wiper material can be provided in pre-wetted formats for reduced VOC emissions, increased convenience, increased productivity, improved solvent control and cleaning protocol repeatability and reduced costs. Pre-wetted options include WFI .

Validated Sterile

ANSI/AAMI/ISO procedures have been followed for verification of irradiation dose with audits being performed to monitor bioburden levels and dose. The Sterility Assurance Level (SAL) or the probability of a viable organism being present on a product unit after sterilization is a SAL of 10⁻⁶ (which is the chance of 1 non-sterile unit in 1,000,000).

Other Sterile wipers

- Gamma Wipe SuperSorb®
- Gamma Wipe® 300
- Gamma Wipe® 120

www.berkshire.com

Contact: Tel 1 800 242 7000 / 1 413 528 2602
info@berkshire.com

America Tel 1 413 528 2602

info@berkshire.com

Europe Tel + 44 1953 562800

enquiries@berkshire.uk.com

SE Asia Tel 65 6252 4313

enquiries@berkshire.com.sg

Japan Tel 81 3 4530 9883

master@berkshire.co.jp

Technical Data:

| Attribute | | Units | Value | Test Method |
|-----------------------------|------------------|-----------------------------------|---------------|--|
| Basis Weight | | g/m ² | 68.0 | TAPPI T-410 |
| Caliper | | µm | 264 | TAPPI T-411 |
| Fibers | ≥100µm | fibers/cm ² | 160 | IEST-RP-CC004.4 Sec 7.1.3/Sec 7.2.2 modified |
| Particles | ≥0.5µm | x10 ⁹ /cm ² | 10 | IEST-RP-CC004.4 Sec 7.1.3/Sec 7.2.1 modified |
| Sorbency | Capacity | mL/m ² | 320 | IEST-RP-CC004.4 Sec 9.1 / Sec 9.2 modified |
| | Efficiency | mL/g | 4.7 | |
| | Rate | seconds | 2 | |
| Non-Volatile Residue | DI Water | g/m ² | 0.028 | IEST-RP-CC004.4 Sec 8.1.2 |
| | IPA | g/m ² | 0.0038 | |
| Ions | Na ⁺ | ppm | 62 | IEST-RP-CC004.4 Sec 8.2.2 |
| | K ⁺ | ppm | 5.9 | |
| | Ca ⁺⁺ | ppm | 22 | |
| | Mg ⁺⁺ | ppm | 5.0 | |
| | Cl ⁻ | ppm | 31 | |

Notes:

- Technical data represented in this table are typical values at the time of publication. These should not be used as product specifications.
- Due to differences in test methods applied and equipment utilized by different wiper manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions
- Third party testing can be performed upon request

Order Information:

| Product | Number | Size | Shts/pk | Pks/Inner bag | Inner bags/cs | Pks/cs | Style |
|----------------|------------|-----------------|---------|---------------|---------------|--------|------------|
| Gamma Wipe® 67 | GW67ST0980 | 9x9"(23x23cm) | 20 | 5 | 16 | 80 | Header Bag |
| Gamma Wipe® 67 | GW67ST1242 | 12x12"(30x30cm) | 25 | 7 | 6 | 42 | Header Bag |
| Gamma Wipe® 67 | GW67ST25 | 12x12"(30x30cm) | 20 | 5 | 5 | 25 | Header Bag |
| Gamma Wipe® 67 | GW67ST60 | 12x12"(30x30cm) | 5 | 10 | 6 | 60 | Header Bag |

Other Berkshire Products



Wipers



Glove Liners



Mop Systems



Documentation Systems



Face Masks



Swabs